

Baby Bath Thermometer

A great gift for the new mother

Safe Accurate Unbreakable Durable Easy to use

Enhances the understanding of the importance of the temperature of the baby bath. Comes with easy instructions on the rear and photo showing how to hold the baby safely.

The Baby Bath Thermometer is an ideal promotional item:

- 1. Available space for printing.
- 2. Long product life as the child when they get older, will insist of using it every time they have a bath.

Baby Bath Thermometer 50° 0 48° A M 46° **Baby Bath** G **Thermometer** 440 Instructions for use 42° 1. Stir bath water to R get even temperature 40° Immerse thermometer fully into the water and keep it there for 15 seconds. 38° IDEAL Temperature must be read 36° while the thermometer is still under water. 340 4. Read the brightest number. If two TOO numbers are equally bright, the 32° temperature is in-between. COLD Note: If the thermometer stays 30° black, the water is either too hot or too cold. Care of thermometer Do not leave in direct sunlight. Wipe dry after use.

3. Flat, lightweight and easy to mail in a standard letter.

Instructions: Stir water to get even temperature. Put thermometer fully into water and keep it there for 15 seconds. Temperature must be read while thermometer is still under water.

Reading: The liquid crystals light up when their respective temperatures are reached. Read the brightest number. If two numbers are equally bright, the temperature is in between.

Care of the thermometer:

Do not leave in direct sunshine. Wipe dry after use

Specifications:

Temperature range:	+30 to +50°C	Dimensions:	146 x 66mm
Resolution:	2°C Increments	Weight:	3.8 gm
Accuracy:	±0.5°C		

Temperature Technology Tel: 08-8297 7077 email: sales@t-tec.com.au 105 Anzac Highway Fax: 08-8371 0558 Web: www.t-tec.com.au Ashford SA 5035